

## Competent Cell Preparation

1. **5 ml** of **LB Medium** where inoculated with *B. Firmus* and *B. Subtilis* respectively
2. The cells where incubated **over night** at 37°C and 200 rpm
3. **3.125 ml** of the ON cultures where transferred into **50 ml** of **growth medium** (diluted 16-Fold)
4. The culture was incubated at 37°C and 200 rpm until the OD<sub>600</sub> reached 1.0  
OD Measurement *B. Firmus*: 0.404 – 0.456 – 0.848 – 1.286 (Avarage of 5 values)  
OD Measurement *B. Subtilis*: 0.440 – 0.490 – 0.816 – 1.310
5. The culture was chilled by transferring the tube on ice and keeping it there for **10 minutes**
6. The culture was centrifuged at 3202 rcf for **30 minutes** at 4°C
7. The pellet was suspended in **15 ml** of **electroporation (washing-) buffer** precooled to 4°C
8. Centrifugation and suspension was repeated **3 times**
9. The pellet of the last centrifugation was suspended in **1.25 ml electroporation buffer**  
OD Measurement *B. Firmus*: 40.8  
OD Measurement *B. Subtilis*: 36.34
10. The OD was set to 100 by centrifugation (5000 rcf / **5 min**) and resuspending in the calculated volume of **electroporation buffer supernatant** to achieve an OD of 100
11. The cells where stored at -80°C to perform electroporation the next day

## Electroporation

1. **60 µl** of thawed (on ice) cells where transferred to precooled electroporation cuvettes with 1 mm gap
2. **100 ng precooled DNA** was added to the electroporation cuvette
3. The mixture inside the gap was stirred with the pipette and gently knocking of the cuvette on the surface of the bench covered with tissue
4. The mixture was incubated for **15 minutes** on ice
5. The mixture was shocked with a single pulse (**3.5 – 5 ms**)
6. The cells where diluted with **1 ml** of **recovery medium** immediately after the shock
7. The recovering cells where incubated for **3 hours** at 37°C and 200 rpm in 15 ml falcon tubes positioned in an angle to increase the surface of the mixture to the surrounding oxygen
8. The cells where spread on selective medium and non-selective medium

## Samples

<u>Organism</u>	<u>Plasmid</u>	<u>Antibiotic</u>
1. Bacillus Subtillis	pAD43-25	Chloramphenicol
2. Bacillus Subtillis	pHT01	Chloramphenicol
3. Bacillus Subtillis	pBSMul1	Kanamycine
4. Bacillus Subtillis	Negative Control	Chloramphenicol
5. Bacillus Firmus	pAD43-25	Chloramphenicol
6. Bacillus Firmus	pHT01	Chloramphenicol
7. Bacillus Firmus	pBSMul1	Kanamycine
8. Bacillus Firmus	Negative Control	Kanamycine

## Troubleshooting and later changed parameters

- Ratio of cell concentration (OD 100) to the mass of plasmid (100 ng)?
- Centrifugation of recovered cells and resuspending in volume-reduced supernatant after 3 h incubation to increase the chance of having transformed cells on one regular sized selective LB agarplate. The total volume of the recovered mixture is 60 µl cells (+ 2-5 µl DNA) + 1 ml recovery medium = 1.062 – 1.065 ml. The protocol advises to use a maximum of 200 µl per agarplate, requiring 5 agarplates per sample.
- Calculations correct ? (see Page 5)
- Applied Voltage, booth 1800 KV and 2100 KV where tried with electroporation cuvettes with 1 mm gap with a time constant of at least 5 ms
- Time of Incubation after electroporation: 3h and 4h
- Cells where concentrated before electroporation to OD 120 in a later performed electroporation
- Weight of used Plasmid: 200 ng of plasmid DNA instead of 100 ng was mixed with the electrocompetent cells
- Presence of organism of interest *B.firmus* on the currently used plates and the positive controls after electroporation was verified via colony PCR
- The protocol worked for *E.coli*: Cells where able to grow on ampicillin selective medium where cells without plasmid where not growing. Gfp expression was verified by fluorescence microscopy and plasmid DNA was harvested via plasmid preparation after preparation of an over Night culture of transformed *E.coli* cells in ampicillin selective liquid medium

## Media

### Growth Medium (55 ml + Blank)

LB: 1% Tryptone, 0,5% yeast extract, 0,5% NaCl  
Sorbitol: 0,5 mol / L  
Glycine Betaine: 7,5 % (W/V)

### Electroporation Buffer (0,5 ml)

Sorbitol: 0,5 mol / L  
Mannitol: 0,5 mol / L  
Glycine Betaine: 7,5 % (W/V)  
Glycerol: 10 % (W/V)

### Electroporation (Washing-) Buffer (150 ml)

Sorbitol: 0,5 mol / L  
Mannitol: 0,5 mol / L  
Glycerol: 10 % (W/V)

### Recovery Medium (1 ml)

LB: 1% Tryptone, 0,5% yeast extract, 0,5% NaCl  
Sorbitol: 0,5 mol / L  
Mannitol: 0,38 mol / L

## Calculations

### Growth Medium 100 ml

Sorbitol:

$$V = 100 \text{ ml} = 0,1 \text{ l}$$

$$c = 0,5 \text{ mol / l}$$

$$M = 182,2 \text{ g / mol}$$

$$n = c * V = 0,5 \text{ mol / l} * 0,1 \text{ l} = 0,05 \text{ mol}$$

$$m = M * n = 182,2 \text{ g / mol} * 0,05 \text{ mol} = \underline{9,11 \text{ g}}$$

Glycine Betaine: (7,5 % (W/V)):

$$m1 / V1 = m2 / V2$$

$$m2 = m1 * V2 / V1 = 7,5 \text{ g} * 100 \text{ l} / 100 \text{ l} = \underline{7,5 \text{ g}}$$

### Electroporationbuffer 10 ml

Sorbitol:

$$V = 10 \text{ ml} = 0.01 \text{ l}$$

$$c = 0,5 \text{ mol / l}$$

$$M = 182,2 \text{ g / mol}$$

$$n = c * V = 0,5 \text{ mol / l} * 0,01 \text{ l} = 0,005 \text{ mol}$$

$$m = M * n = 182,2 \text{ g / mol} * 0,005 \text{ mol} = \underline{0.91 \text{ g}}$$

Mannitol:

$$V = 10 \text{ ml} = 0.01 \text{ l}$$

$$c = 0,5 \text{ mol / l}$$

$$M = 182,2 \text{ g / mol}$$

$$n = c * V = 0,5 \text{ mol / l} * 0,01 \text{ l} = 0,005 \text{ mol}$$

$$m = M * n = 182,2 \text{ g / mol} * 0,005 \text{ mol} = \underline{0.91 \text{ g}}$$

Glycine Betaine: (7,5 % (W/V)):

$$m1 / V1 = m2 / V2$$

$$m2 = m1 * V2 / V1 = 7,5 \text{ g} * 10 \text{ ml} / 100 = \underline{0.75 \text{ g}}$$

Glycerol: 10 % (W/V):

$$c1 * V1 = c2 * V2$$

$$10 \% * 10 \text{ ml} = 87 \% * x \text{ ml}$$

$$x \text{ ml} = 10 \% * 10 \text{ ml} / 87 \% = \underline{1.15 \text{ ml}}$$

$$10 \text{ ml} = 1.15 \text{ ml (Glycerol)} + 8.85 \text{ ml (ddH}_2\text{O)}$$

## Electroporation (Washing-) Buffer 150 ml

### Sorbitol:

$$V = 150 \text{ ml} = 0.15 \text{ l}$$

$$c = 0,5 \text{ mol / l}$$

$$M = 182,2 \text{ g / mol}$$

$$n = c * V = 0,5 \text{ mol / l} * 0,15 \text{ l} = 0,075 \text{ mol}$$

$$m = M * n = 182,2 \text{ g / mol} * 0,075 \text{ mol} = \underline{13.66 \text{ g}}$$

### Mannitol:

$$V = 150 \text{ ml} = 0.15 \text{ l}$$

$$c = 0,5 \text{ mol / l}$$

$$M = 182,2 \text{ g / mol}$$

$$n = c * V = 0,5 \text{ mol / l} * 0,15 \text{ l} = 0,075 \text{ mol}$$

$$m = M * n = 182,2 \text{ g / mol} * 0,075 \text{ mol} = \underline{13.66 \text{ g}}$$

### Glycerol: 10 % (W/V):

$$c_1 * V_1 = c_2 * V_2$$

$$10 \% * 150 \text{ ml} = 87 \% * x \text{ ml}$$

$$x \text{ ml} = 10 \% * 150 \text{ ml} / 87 \% = 17.24 \text{ ml}$$

$$150 \text{ ml} = 17.24 \text{ ml (Glycerol)} + 132.76 \text{ ml (ddH}_2\text{O)}$$

---

## Recovery Medium 50 ml

### Sorbitol:

$$V = 50 \text{ ml} = 0,05 \text{ l}$$

$$c = 0,5 \text{ mol / l}$$

$$M = 182,2 \text{ g / mol}$$

$$n = c * V = 0,5 \text{ mol / l} * 0,05 \text{ l} = 0,025 \text{ mol}$$

$$m = M * n = 182,2 \text{ g / mol} * 0,025 \text{ mol} = \underline{4,555 \text{ g}}$$

### Mannitol:

$$V = 50 \text{ ml} = 0,05 \text{ l}$$

$$c = 0,38 \text{ mol / l}$$

$$M = 182,2 \text{ g / mol}$$

$$n = c * V = 0,38 \text{ mol / l} * 0,05 \text{ l} = 0,019 \text{ mol}$$

$$m = M * n = 182,2 \text{ g / mol} * 0,019 \text{ mol} = \underline{3,4618 \text{ g}}$$

## Antibiotics

### Chloramphenicol

Stock: 25 mg / ml

Work: 10 µg / ml

V = 500 mL (= 20 Plates)

$$m_{\text{work}}: m_1 / V_1 = m_2 / V_2$$

$$\leftrightarrow 10 \mu\text{g} / 1 \text{ ml} = X / 500 \text{ ml}$$

$$\leftrightarrow X = 10 \mu\text{g} * 500 \text{ ml} / 1 \text{ ml} = \underline{5 \text{ mg}} \text{ (CM required to have 500 ml of a } 10 \mu\text{g} / \text{ml solution)}$$

$$V_{\text{stock}}: m_1 / V_1 = m_2 / V_2$$

$$\leftrightarrow 25 \text{ mg} / 1 \text{ ml} = \underline{5 \text{ mg}} / X$$

$$\leftrightarrow X = \underline{5 \text{ mg}} * 1 \text{ ml} / 25 \text{ mg} = \underline{0,2 \text{ ml}} \text{ (of } 25 \text{ mg} / \text{ml CM stock solution is required for 500 ml of a } 10 \mu\text{g} / \text{ml work solution)}$$

### Kanamycine

Stock: 50 mg / ml

Work: 50 µg / ml

V = 500 mL (= 20 Plates)

$$m_{\text{work}}: m_1 / V_1 = m_2 / V_2$$

$$\leftrightarrow 50 \mu\text{g} / 1 \text{ ml} = X / 500 \text{ ml}$$

$$\leftrightarrow X = 50 \mu\text{g} * 500 \text{ ml} / 1 \text{ ml} = \underline{25 \text{ mg}} \text{ (KM required to have 500 ml of a } 50 \mu\text{g} / \text{ml solution)}$$


$$V_{\text{stock}}: m_1 / V_1 = m_2 / V_2$$

$$\leftrightarrow 50 \text{ mg} / 1 \text{ ml} = \underline{25 \text{ mg}} / X$$

$$\leftrightarrow X = \underline{25 \text{ mg}} * 1 \text{ ml} / 50 \text{ mg} = \underline{0,5 \text{ ml}} \text{ (of } 50 \text{ mg} / \text{ml KM stock solution is required for 500 ml of a } 50 \mu\text{g} / \text{ml work solution)}$$


Plasmids:

PBSMuL




P<sub>HpaII</sub>, constitutive Gram-positive promoter;  
P<sub>59</sub>, second constitutive gram-positive promoter only present in pBSMuL2;  
MCS, multiple cloning site;  
ori E. coli, gene required for replication in *E. coli*;  
amp<sub>r</sub>,  $\beta$ -lactamase gene conferring ampicillin resistance in *E. coli*;  
repB, gene required for replication in *B. subtilis*;  
kan<sub>r</sub>, kanamycin resistance gene

pHT01



The pHT01 is a high-level intracellular expression vector for recombinant protein production with *B. subtilis*. The vector is based on the strong  $\sigma^A$ -dependent promoter preceding the *groES-groEL* operon of *B. subtilis* which has been converted into an efficiently controllable (IPTG-inducible) promoter by addition of the *lac* operator. pHT01 is an *E. coli*/*B. subtilis* shuttle vector, that provides ampicillin resistance to *E. coli* and chloramphenicol resistance to *B. subtilis*.

## PAD43-25



- Features:**
- gfpmut3a* promoter-less gene encoding a variant of Green Fluorescent Protein from plasmid pFPV25 Valdivia, R. H. and S. Falkow. 1997. *Science* **277**:2007-2011
  - rep* Replication initiation protein from cryptic rolling circle plasmid pTA1060 (GenBank U32380) from *Bacillus subtilis* "natto"
  - cat* encodes chloramphenicol acetyl transferase; selectable in either *E. coli* or *B. subtilis* (chloramphenicol 5 µg/ml)
  - bla* encodes β-lactamase; selectable in *E. coli* only (ampicillin 100 µg/ml)
  - glyA* last 1010 bp of the *Bacillus cereus* structural gene for glycine/serine hydroxymethyltransferase
  - P<sub>upp</sub>* constitutive promoter from the *Bacillus cereus upp* (uracil phosphoribosyltransferase) gene

**Description:** pAD43-25 is a shuttle vector, replicating in *E. coli* from the pBR322 origin and in *Bacillus* from the pTA1060 origin. A chromosomal fragment from *Bacillus cereus* UW85, containing the first 21 bp of the *upp* gene and all of its upstream regulatory regions, allows for high-level constitutive expression of a Green Fluorescent Protein variant.

**Construction:** pAD43-25 was constructed by ligating chromosomal DNA *Sau3A* fragments from UW85 into the pAD123 *Bam*HI site and screening for clones that express high levels of GFP in *Bacillus cereus* during vegetative growth.

**Use:** The placement of the *B. cereus* UW85 *upp* promoter upstream of *gfpmut3a* allows for constitutive expression of a mutant GFP that has been optimized for use in fluorescence-activated cell sorting, with an optimal excitation wavelength of 498 nm. This shuttle vector should replicate in a wide variety of Gram-positive organisms along with *E. coli*. Plasmid pAD43-25 should serve as a useful marker for quantifying or sorting cells in a wide variety of applications.